TECHNICAL SPECIFICATIONS

SOPIX²

Size 1

External dimensions	25 x 39mm
Active surface area	600mm2 (20 x 30mm)
Number of pixels	1.50million

SOPIX inside system

Technology	CMOS + scintillator+ optic fiber
	20µm x 20µm
	25lp/mm Real
	>12lp/mm Supplied
	Sopro Imaging TWAIN
module	Yes

SOPIX² USB connection

Connection	USB 2.0 Total
cable length	3 70m

External dimensions	31 x 42mm
Active surface area	884mm2 (26 x 34mm)
Number of pixels	2.21millions

SOPIX²

Technology	CMOS + scintillator + optic fiber
Pixel size	20μm x 20μm
Theoretical resolution	25lp/mm
Real resolution	>18lp/mm
Supplied imaging softwa	reSopro Imaging
TWAIN module	Yes

SOPIX inside USB connection

Connection	USB 2.0
Sensor cable length	0.70m

PSPIX²

System

Resolution	20 lp/mm
Scan Time (fast mode)	1,6s - 2,7s
Scan Time (high definition mode)	2,1s - 3,6s
Connection	Ethernet RJ-45
Dimensions L. 154 x [D. 204 x H. 193 mm
Weight	2,6 kg
Operating voltage100	

Imaging Plates

Dimensions IP Size 0	22 x 35 mm
Dimensions IP Size 1	
Dimensions IP Size 2	31 x 41 mm
Dimensions IP Size 3	27 x 54 mm
Dimensions IP Size 4 (3 x IP Size 3)	69 x 54 mm
(

Windows® minimum configuration required Operating system......Windows 10

MAC® minimum configuration required

Hard disk......1 TB Graphic card...... Chipset Nvidia® or ATI® 2 GB unshared memory compatible DirectX 9 or more Screen resolution.....1280 x 1024 or more Ethernet board......1 Gbps

Operating system...... OS X El Capitan

COMPUTER CONFIGURATION

Windows® minimum configuration required

Operating system	Windows / SP1
Processor	Core 2 Duo - 3GHz
RAM	2 GB
Hard disk	250 GB
Graphic card	512 MB RAM unshared memory
	compatible DirectX 9
Screen resolution	1280 x 1024
Ethernet board	100 Mbps - 1 Gbps
MAC® minimum	configuration required
Computer	MacBook® Pro 13.3" or iMac® 21.5"
	OS X Mavericks
	Intol® Coro 2 Duo

For Yosemite and El Capitan operating systems, a Mac computer from 2013 or later is required.

Note: In the case of SOPIX inside and SOPIX2 inside, the IEC 60601-2-65 norm requires for each X-Ray intraoral system with an onboard digital sensor to use a

Note: The data transfer from the intraoral system X-Mind unity to SOPRO® Imaging is not available on SOPRO® Imaging Mac version yet.

The medical devices for dental care SOPIX® Series are of class IIa and manufactured by SOPRO®, notified body LNE/GMED, X-Mind unity is of class IIb and manufactured by DE GOTZEN, notified body DNV - CE 0434. These medical devices are not refunded by health insurance organizations. Read carefully the instructions on the labelling before use.


The medical device for dental care PSPIX® is of class IIa and manufactured by SOPRO® notified body LNE/GMED. This medical device is not refunded by health insurance organizations. Read carefully the instructions on the labelling before use.

PSPIX® and SOPRO® are registered trademarks of SOPRO. "All other trademarks cited herein are the property of their respective owners"

124 Gaither Drive | Suite 140 Mount Laurel | NJ | 08054 | USA Fax. 1-856-222-4726 www.acteongroup.com info@acteonusa.com

MORE INVENTIVE LESS INVASIVE





The cordless imaging plate scanner



SOPIX

A perfect image every time with minimum exposure to radiation










STRIKING CONTRAST FOR A MORE RELIABLE DIAGNOSIS









MORE INVENTIVE

Better differentiation of dental tissue

SOPIX® sensors surpass the limits of radiological examinations by offering greater differentiation of dental tissue.

This technological achievement is based on the use of **broad spectrum optical** microfibers for the guided transmission of photon emissions in order to provide **highly contrasted images.**





LESS INVASIVE

A more reliable diagnosis

The different tooth anatomic structures such as the bone, roots, pulp... are highlighted with **extreme precision** in the image.

Your diagnosis is **faster** and **more accurate**!





A perfect image every time with minimum exposure to radiation



THE PERFECT FIT FOR YOUR CLINICAL PRACTICE




Endodontics





Bitewing



Pediatric



Periapical



Implantology



Cariology



SOPIX sensors provide accurate images and striking contrast to ensure a reliable diagnosis.



Two sizes are available depending on patient morphology and clinical applications.







A QUALITY IMAGE

VIA AN INTERFACE THAT IS

The ACTEON® Imaging Suite software offers intuitive navigation with the mouse and advanced functionality. It alone lets you manage all of your images, from scanning to viewing images from all ACTEON® imaging devices (CBCT, Panoramic, intraoral digital X-ray system, intraoral camera, etc.) and much more.





A QUALITY IMAGE EVERY TIME WITH MINIMAL EXPOSURE TO RADIATION


CUTTING EDGE TECHNOLOGY SOPIX

Available in all **SOPIX** series sensors, patented \bigcirc technology (**Automatic Control Exposure**) analyzes in real-time, the amount of X-rays accumulated by the sensor.

It automatically freezes the image acquisition as soon as the sensor receives the radiation required to produce the perfect image.

Eliminate the risk of over exposing the image!

Combined with the **X-Mind unity** intraoral X-ray generator, **SOPIX**® **inside** with **ACE**® **technology limits the emission of x-rays** during the acquisition to the necessary amount for the patient's morphology. It uses the minimum dose required to provide a high-quality image.



Laurent **GUILHAUMON** R&D Project Manager, HW/Embedded SW Systems "ACE is the combination of advanced sensor technology, digital power electronics and the know-how of two diagnostic imaging divisions. The synergy between La Ciotat (FRANCE) and Milan (ITALY) R&D teams gave birth to an innovative concept focused on patients, with outstanding image quality. FOR A SAFER PROCESS

With **SOPIX** digital sensors and its patented **ACE**® **technology**, you acquire **successful X-rays every time**, meaning reliable and accurate diagnosis. You will save time by avoiding the need for retakes.

Using **X-Mind unity** intraoral X-ray generator with **SOPIX inside**, the patients receive the minimum dose required for their dental morphology. You protect your patients and your staff from any unnecessary radiation exposure.

ACTEON Imaging Suite records the **X-Mind unity** settings as well as the effective dose received by the patient. This ensures permanent traceability for every patient.

PATIENT AND STAFF

PROTECTION





ACTEON Imaging Suite, always one step ahead

ACTEON Imaging Suite systematically records the **X-Mind unity** settings as well as the effective dose received by the patient for each acquisition.

This ensures **permanent traceability** for every patient



EXCLUSIVE TRACEABILITY



Outstanding working comfort

Through direct integration of **SOPIX inside** into **X-Mind unity**, connecting cables are hidden inside the X-ray unit.


The holder places the **sensor safely, within reach,** to prevent it from falling onto the floor.

Your working environment is **more ergonomic and productive.**




When **SOPIX inside** has received enough energy to provide an **exceptional image quality**, it tells

The patient only receives the necessary dose adapted for their dental morphology, which **protects them from unnecessary overexposure.**



the X-Mind unity to stop the X-ray emission.

STRIKING CONTRAST PSP!X2 FOR AN ACCURATE DIAGNOSIS




MORE INVENTIVE

Better differentiation of dental tissue

PSPIX² surpasses the limits of radiological examinations by offering greater differentiation of dental tissue.

This technological achievement is called FIBER2PIXEL®

WITH FIBER - FIBER2PIXEL®







Differentiation of dental tissue

FIBER2PIXEL® technology is based on the use of broad spectrum optical microfibers for the guided transmission of photon emissions in order to provide highly contrasted images.

LESS INVASIVE

A more reliable diagnosis

The different tooth anatomic structures, such as the bone, roots, pulp... are highlighted with **extreme precision** on the image.

Your diagnosis is **faster** and **more accurate!**









A SCANNER FOR MYSELF

ACTEON INTRODUCES THE FIRST CORDLESS DIGITAL IMAGING SCANNER

-The PSPIX² is-

SMALL – compact size fits in any operatory

INTUITIVE – easy to learn and easier to use

ELEGANT – a perfect solution for digital imaging

AFFORDABLE - every chairside can be equipped for imaging





I WANT A SIMPLER AND MORE INTUITIVE SCANNER

PSPIX² IS SO INTUITIVE, SCANNING AN IMAGE HAS NEVER BEEN EASIER

Created with cutting edge technology and design, the PSPIX² is the smallest and most compact PSP scanner on the market.

The large color touchscreen provides simple instructions for quick and easy use.


With clean lines and an elegant design, the **PSPIX²** fits perfectly into any practice.

Workflow has never been so smooth and efficient, through its outstanding intellectual ability. A real evolution in your working practice!



PSPIX² features optional removable parts for easy disinfection. The autoclavable parts further provide a high level of protection


PSPIX² provides accurate, sharp and contrasted images to ensure reliable clinical diagnosis












PSPIX² is delivered with **ACTEON Imaging Suite** software, which is easy to integrate and used by thousands

ACTEON Imaging Suite is compatible with MAC® and Windows® operating systems



MAC COMPATIBLE





I CAN SCAN EVERY PATIENT

PSPIX² EXCEEDS YOUR TREATMENT EXPECTATIONS

PSPIX²



With increased flexibility, the wireless **ACTEON** imaging plates provide greater comfort for all your patients.

- Ideal for both children and adults
- Imaging plates are available in sizes: 0, 1, 2, 3
- Endorsed by pediatric dentists worldwide







CAN SHARE IT

PSPIX² IS ALSO DESIGNED FOR MULTIUSER ENVIRONMENTS

LET THE LIGHT GUIDE YOU...

Check the status of **PSPIX²**at a quick glance!



Blue: Available



Purple: Scanning



Yellow: Occupied

JUST CLICK AND SCAN!

With award winning **ACTEON Imaging Suite** software the **PSPIX**² can be shared across your operatories. Simply reserve the scanner or select your workstation on the touchscreen, insert your plate and let the **PSPIX**² do the rest... Images scan in seconds and are immediately available.



FULLY AUTOMATED PROCESS

- AUTO-ACCESS door only opens when a plate is detected
- AUTO-DETECT plate size is automatically detected
- AUTO-SCAN PSPIX² scans and optimizes images before saving them into a patient record
- AUTO-EJECT plates are automatically erased and ejected
- **ECO-MODE** enters standby-mode to save power when not in use









